

RNA extraction with Trizol

For Whole Tissues:

- 1 – Isolate tissue and snap freeze in liquid nitrogen
- 2 – Transfer tissue to a 14-ml culture tube and add 2ml of Trizol
- 3 – Homogenize completely
- 4 – Transfer 1ml each into two fresh 1.5-ml Epp. tubes
- 5 – Incubate 5 min @ RT
(Recommended: Freeze in Trizol at least 1 hr @ -80°C after incubating @ RT)
*Potential indefinite stopping point @ -80°C

For Sorted Cells:

- 1 – Sort cells into 1.5 ml epp. tubes coated in Fetal Calf Serum (FCS) or in 500µl media with 2% FCS
(Recommended: For ≤ 100,000 cells sort directly into 500µl Trizol and proceed to step 5)
- 2 – Centrifuge 30 sec @ full speed, 4°C (make sure to see a pellet)
- 3 – Remove supernatant from each sample with a pipette tip
Leave behind a small amount of liquid so as not to disturb the pellet
- 4 – Add 500 µl Trizol to each sample (mix well by pipetting)
- 5 – Incubate 5 min @ RT
(Recommended: Freeze in Trizol at least 1 hr @ -80°C after incubating @ RT)
* Potential indefinite stopping point @ -80°C

Thaw Trizol samples at RT

- 6 – Add the appropriate amount of Chloroform to each sample (work in hood, do not pipette)
 - For 500µl Trizol, Add 100µl Chloroform
 - For 1ml Trizol, Add 200µl Chloroform
- 7 – Shake vigorously for approx 15 sec
- 8 – **Incubate 2-3 min @ RT**
- 9 – Centrifuge **15 min @ 12,000 g, 4°C**
- 10 – Transfer aqueous phase CAREFULLY to a fresh tube¹ (work on ice from this point on)
 - For 500µl Trizol, Transfer ~250µl
 - For 1ml Trizol, Transfer ~500µl
- 11 – Add 1µl Glycoblue to each sample (flick well but don't pipette to mix, quick spin, on ice)
- 12 – Add the appropriate amount of Isopropanol to each sample (flick well to mix, do not pipette)
 - For 500µl Trizol, Add 250µl Isopropanol
 - For 1ml Trizol, Add 500µl Isopropanol
- 13 – **Incubate 1+ hr @ -80°C**
*Potential indefinite stopping point @ -80°C
- 14 – Centrifuge **20 min @ 12,000 g, 4°C**
Remove supernatant with a pipette tip
- 15 – Wash with the appropriate amount of 75% cold ethanol (flick to loosen pellet)
 - For 500µl Trizol, Wash with 500µl 75% cold ethanol
 - For 1ml Trizol, Wash with 1ml 75% cold ethanol
- 16 – Centrifuge **15 min @ 7,400 g, 4°C**
Pipette off as much ethanol as possible, air dry 1 hr (pellet will change color)
- 17 – Resuspend each RNA pellet in 12.5 - 15 µl Nuclease-free water
- 18 – Incubate @ RT for **2-3 min**, then flick well to mix and place on ice
- 19 – Quantify total RNA using a Nanodrop
*Store total RNA indefinitely @ -80°C

¹ Leave behind a small amount of clear aqueous phase; do NOT pick up any pink phenol-chloroform phase; use pipette tips with a larger hole to prevent this from happening.