

Figure S1


A

Immunological Genome
Project
ImmGen.org

Edy Kim
Brenner Lab
Brigham & Women's Hospital

NKT.Sp.LPS.3hr: Splenic iNKT, 3h post-LPS challenge

Pre-sort



Post-1st sort

Tissue Preparation : Mechanical disruption of spleen


Time from Mouse to Trizol: 2.5 hours

Antibody (Clone) - Fluorochrome

TCRb (H57-597)—V450

CD8 (53-6.7), CD19 (1D3), Gr1 (8C5), Ter119 (Ter119) -- APC-Cy7


CD1d tetramer loaded with PBS-57 -- APC
(provided by NIH tetramer facility)




94.6 pure
after 1st
sort, >>95%
pure on 2nd
sort

Treg.4.25hi.Sp: Splenic regulatory T cells

Pre-sort



2nd sorting



Tissue Preparation:

100um mesh/syringe plunger teasing, ACK lysis,


Time from Mouse to the 2nd sorting:

8am to 10am

Antibody (Clone) - Fluorochrome:

CD4 (RM4-5) – APC, CD8a (53-6.7) – A700, TCRb (H57-597) – eF450,
CD44(IM7) – FITC, CD62L(Mel-14) – PE-Cy7, CD25(PC61) -- PE
dump = CD19(1D3), Ter119(Ter-119) , Gr1(RB6-8C5) – APC-eF780
NK1.1(PK136), F4/80(BM8) -- APC-Cy7

Treg.4.FP3+.Nrpl0.Co: Colonic regulatory T cells



Tissue Preparation:

removing fat, stripping IEL layer and mucus by EDTA+DTT,
digesting by dispase+collagenase


Time from Mouse to the 2nd sorting:

7am to 12am


Antibody (Clone) - Fluorochrome:

CD4 (RM4-5) – PE, TCRb (H57-597) – PE-Cy7, Foxp3-IRES-GFP,
Nrpl0 – APC
dump = CD19(1D3) – PE-TxRed

B.Fem.SP: Splenic B cell (Female)



Post-sort



Tissue Preparation: Mechanical disruption of spleen

Time from Mouse to TCL Buffer: 10:30am-12:00pm


Antibody (Clone) - Fluorochrome:

TCRb (H57-597)—Pacific Blue


CD19 (1D3)—PE-Cy7

IgM (eB121-15F9)—FITC

B.Sp: Splenic B cells



2nd sorting



Tissue Preparation:

100um mesh/syringe plunger teasing, ACK lysis

Time from Mouse to the 2nd sorting:

9:30am to 11am


Antibody (Clone) - Fluorochrome:

TCRb(H57-597)--PB

CD19(1D3)-- PE-Cy7

IgM(eB121-15F9)--FITC

Ep.MECHi.Th: Thymic medullary epithelial cells



Post-sort

Tissue Preparation:

Collagenase digestion, Collagenase/Dispase digestion, CD45 MACS depletion (LS column)

Time from Mouse to the 2nd sorting:

8am to 12am

Antibody (Clone) - Fluorochrome:


CD45 (30-F11) -- PE-CF594

Ly51 (6C3) -- PE


MHCII(M5/114) -- FITC

EpCAM(G8.8) -- APC

MF.Fem.PC: Peritoneal cavity macrophages (female)



2nd sorting



Tissue Preparation:

Peritoneal lavage with Phenol-red free DMEM, ACK lysis

Time from Mouse to the 2nd sorting:

8am to 10am

Antibody (Clone) - Fluorochrome:


F4/80 (BM8) -- PE-Cy7

ICAM2 (3C4 (mIC2/4)) -- FITC


MF.PC: Peritoneal cavity macrophages

1) Single Cells


Pre-sort




2) Live Cells



3) MF: F4/80+ ICAM2+



Post-sort



Tissue Preparation: Peritoneal Cavity Lavage


Time from Mouse to TCL Buffer: 8:30am-11:00am

Antibody (Clone) - Fluorochrome:


ICAM2 (mIC2/4)—FITC

F4/80 (BM8)—PE

preT.DN1.Th: Pre-TCR DN1 thymocytes



2nd sorting



Tissue Preparation:

100um mesh/syringe plunger teasing, ACK lysis,
MACS depletion with APC-beads and LS column

Time from Mouse to the 2nd sorting:

8am to 11am

Antibody (Clone) - Fluorochrome:

CD25 (PC61) – PE, CD44 (IM7) – FITC


CD117 (2B8) – PE-Cy7

Lin = CD4(RM4-5), CD8a(53-6.7), NK1.1(PK136) – APC


CD11b(M1/70), TCRb(H57-597), TCRgd(GL3) – PerCP-Cy5.5

CD11c(N418), CD19(1D3), Ter119(Ter-119), Gr1(RB6-8C5) – APC-Cy7

preT.DN2a.Th: Pre-TCR DN2a thymocytes



2nd sorting



Tissue Preparation:

100um mesh/syringe plunger teasing, ACK lysis,
MACS depletion with APC-beads and LS column

Time from Mouse to the 2nd sorting:

8am to 11am

Antibody (Clone) - Fluorochrome:

CD25 (PC61) – PE, CD44 (IM7) – FITC


CD117 (2B8) – PE-Cy7

Lin = CD4(RM4-5), CD8a(53-6.7), NK1.1(PK136) – APC


CD11b(M1/70), TCRb(H57-597), TCRgd(GL3) – PerCP-Cy5.5

CD11c(N418), CD19(1D3), Ter119(Ter-119), Gr1(RB6-8C5) – APC-Cy7

preT.DN2b.Th: Pre-TCR DN2b thymocytes



2nd sorting



Tissue Preparation:

100um mesh/syringe plunger teasing, ACK lysis,
MACS depletion with APC-beads and LS column

Time from Mouse to the 2nd sorting:

8am to 11am

Antibody (Clone) - Fluorochrome:

CD25 (PC61) – PE, CD44 (IM7) – FITC


CD117 (2B8) – PE-Cy7

Lin = CD4(RM4-5), CD8a(53-6.7), NK1.1(PK136) – APC


CD11b(M1/70), TCRb(H57-597), TCRgd(GL3) – PerCP-Cy5.5

CD11c(N418), CD19(1D3), Ter119(Ter-119), Gr1(RB6-8C5) – APC-Cy7

preT.DN3.Th: Pre-TCR DN3 thymocytes



2nd sorting



Tissue Preparation:

100um mesh/syringe plunger teasing, ACK lysis,
MACS depletion with APC-beads and LS column

Time from Mouse to the 2nd sorting:


8am to 11am

Antibody (Clone) - Fluorochrome:


CD25 (PC61) – PE, CD44 (IM7) – FITC
CD117 (2B8) – PE-Cy7, CD28(E18) -bio-SA -- PETxRed
Lin = CD4(RM4-5), CD8a(53-6.7), NK1.1(PK136) – APC
CD11b(M1/70), TCRb(H57-597), TCRgd(GL3) – PerCP-Cy5.5
CD11c(N418), CD19(1D3), Ter119(Ter-119), Gr1(RB6-8C5) – APC-Cy7

T.4.Nve.Fem.Sp: Splenic Naïve CD4 T cells (female)

Pre-sort



2nd sorting



Tissue Preparation:

100um mesh/syringe plunger teasing, ACK lysis

Time from Mouse to the 2nd sorting:


8am to 10am

Antibody (Clone) - Fluorochrome:


CD4 (RM4-5) – PE, CD8a (53-6.7) – A700, TCRb (H57-597) – eF450,
CD44(1M7) – FITC, CD62L(Mel-14) – PE-Cy7, CD25(PC61) -- APC
dump = CD19(1D3), Ter119(Ter-119) , Gr1(RB6-8C5) – APC-eF780
NK1.1(PK136), F4/80(BM8) -- APC-Cy7

T.4.Nve.Sp: Splenic Naïve CD4 T cells

Pre-sort



2nd sorting



Tissue Preparation:

100um mesh/syringe plunger teasing, ACK lysis,


Time from Mouse to the 2nd sorting:

8am to 10am


Antibody (Clone) - Fluorochrome:

CD4 (RM4-5) – APC, CD8a (53-6.7) – A700, TCRb (H57-597) – eF450,
CD44(IM7) – FITC, CD62L(Mel-14) – PE-Cy7, CD25(PC61) -- PE
dump = CD19(1D3), Ter119(Ter-119) , Gr1(RB6-8C5) – APC-eF780
NK1.1(PK136), F4/80(BM8) -- APC-Cy7

T.4.Sp.aCD3+CD40.18hr: Splenic CD4 T cells, 18hrs post anti-CD3+CD40



2nd sorting



Tissue Preparation:

ip 60ug CD3(KT3) and 100ug CD40(FGK4.5) 18hrs before sac.
100um mesh/syringe plunger teasing, ACK lysis


Time from Mouse to the 2nd sorting:

8am to 10am


Antibody (Clone) - Fluorochrome:

CD4 (RM4-5) – PE, CD8a (53-6.7) – PE-Cy7, TCRb (H57-597) – eF450
dump = CD19(1D3), Ter119(Ter-119), Gr1(RB6-8C5) – APC-eF780
NK1.1(PK136), F4/80(BM8) -- APC-Cy7

T.4.Th: Thymic CD4SP T cells



2nd sorting



Tissue Preparation:

100um mesh/syringe plunger teasing, ACK lysis,

Time from Mouse to the 2nd sorting:

8am to 10am

Antibody (Clone) - Fluorochrome:


CD4 (RM4-5) – PE, CD8a (53-6.7) – APC, TCRb (H57-597) – PB,
CD24(M1/69) -- FITC

dump = CD19(1D3), Ter119(Ter-119) , Gr1(RB6-8C5) – APC-eF780


NK1.1(PK136), F4/80(BM8) -- APC-Cy7

T.8.Nve.Sp: Splenic Naïve CD8 T cells

Pre-sort



2nd sorting



Tissue Preparation:

100um mesh/syringe plunger teasing, ACK lysis,


Time from Mouse to the 2nd sorting:

8am to 10am


Antibody (Clone) - Fluorochrome:

CD4 (RM4-5) – APC, CD8a (53-6.7) – A700, TCRb (H57-597) – eF450,
CD44(IM7) – FITC, CD62L(Mel-14) – PE-Cy7
dump = CD19(1D3), Ter119(Ter-119), Gr1(RB6-8C5) – APC-eF780
NK1.1(PK136), F4/80(BM8) -- APC-Cy7

T.8.Th: Thymic CD8SP T cells



2nd sorting



Tissue Preparation:

100um mesh/syringe plunger teasing, ACK lysis,

Time from Mouse to the 2nd sorting:

8am to 10am


Antibody (Clone) - Fluorochrome:

CD4 (RM4-5) – PE, CD8a (53-6.7) – APC, TCRb (H57-597) – PB,
CD24(M1/69) -- FITC


dump = CD19(1D3), Ter119(Ter-119), Gr1(RB6-8C5) – APC-eF780

NK1.1(PK136), F4/80(BM8) -- APC-Cy7

T.DN4.Th: DN4 thymocytes



2nd sorting



Tissue Preparation:

100um mesh/syringe plunger teasing, ACK lysis,
MACS depletion with APC-beads and LS column


Time from Mouse to the 2nd sorting:

8am to 11am


Antibody (Clone) - Fluorochrome:

CD25 (PC61) – PE, CD44 (IM7) – FITC
CD117 (2B8?) – PE-Cy7, CD28(E18) -bio-SA -- PETxRed
Lin = CD4(RM4-5), CD8a(53-6.7), NK1.1(PK136) – APC
CD11b(M1/70), TCRb(H57-597), TCRgd(GL3) – PerCP-Cy5.5
CD11c(N418), CD19(1D3), Ter119(Ter-119), Gr1(RB6-8C5) – APC-Cy7

T.DP.Th: CD4+ CD8+ double-positive thymocytes



2nd sorting



Tissue Preparation:

100um mesh/syringe plunger teasing, ACK lysis,


Time from Mouse to the 2nd sorting:

8am to 10am


Antibody (Clone) - Fluorochrome:

CD4 (RM4-5) – PE, CD8a (53-6.7) – APC, TCRb (H57-597) – PB,
CD24(M1/69) – FITC, CD69(H1.2F3) – A700
dump = CD19(1D3), Ter119(Ter-119), Gr1(RB6-8C5) – APC-eF780
NK1.1(PK136), F4/80(BM8) -- APC-Cy7

T.ISP.Th: Immature single-positive thymocytes



2nd sorting



Tissue Preparation:

100um mesh/syringe plunger teasing, ACK lysis,

Time from Mouse to the 2nd sorting:

8am to 10am



Antibody (Clone) - Fluorochrome:

CD4 (RM4-5) – PE, CD8a (53-6.7) – APC, TCRb (H57-597) – PB,
CD24(M1/69) -- FITC

dump = CD19(1D3), Ter119(Ter-119) , Gr1(RB6-8C5) – APC-eF780

NK1.1(PK136), F4/80(BM8) -- APC-Cy7

ILC3.NKp46+.SI: Small intestine NKp46+ ILC3



Tissue Preparation: Small intestine cut into small pieces, IEL removed by washing with EDTA, dissociated with Collagenase 4, lymphocyte isolation with Percoll gradient.

Time from Mouse to Trizol: 09:15am-03:00pm



Antibody (Clone) - Fluorochrome:

CD45 (30-F11) – APC-Cy7
CD3e (145-2C11) – FITC

CD19 (eBio1D3) – FITC
Thy1.2 (53-2.1) – APC
TCRb (H57-597) – PE
CD5 (53-7.3) – PE
KLRG1 (2F1) – PE
CCR6 (140706) – BV421
NKp46 (CS96) – biotin – SA-PeCy7

ILC3: 96.8%
NKp46+: 96.3%
purity

ILC3.NKp46-CCR6-.SI: Small intestine NKp46- CCR6- ILC3



Tissue Preparation: Small intestine cut into small pieces, IEL removed by washing with EDTA, dissociated with Collagenase 4, lymphocyte isolation with Percoll gradient.

Time from Mouse to Trizol: 09:15am-02:45pm

Antibody (Clone) - Fluorochrome:

CD45 (30-F11) – APC-Cy7
CD3e (145-2C11) – FITC

CD19 (eBio1D3) – FITC
Thy1.2 (53-2.1) – APC
TCRb (H57-597) – PE
CD5 (53-7.3) – PE
KLRG1 (2F1) – PE
CCR6 (140706) – BV421
NKp46 (CS96) – biotin – SA-PeCy7

ILC3: 98.8%
DN: 72.9%
purity






ILC2.SI: Small intestine ILC2

Gated on:


1) lymphocytes

2) Singlets





3) Singlets



5) CD3-CD19-NK1.1-7-AAD-TCRb- CD5-



6) CD127+Sca1+



Tissue Preparation: Small intestine cut into small pieces, IEL removed by washing with EDTA, dissociated with Collagenase 4, lymphocyte isolation with Percoll gradient.

Time from Mouse to Trizol: 08:45am-01:45pm



Antibody (Clone) - Fluorochrome:

CD45 (30-F11) – APC-Cy7
CD3e (145-2C11) – PerCP-Cy5.5

CD19 (eBio1D3) – PerCP-Cy5.5
NK1.1 (PK136) – PerCP-Cy5.5
TCR α (H57-597) – PE
CD5 (53-7.3) – PE
Sca1 (D7) – PacBlue
CD127 (A7R34) – FITC
KLRG1 (2F1) – APC
ST2 (RMST-2) – biotin – SA-PeCy7

ILC2: 98.4%
ST2⁺: 67.5%
purity

ILC3.CCR6+.SI: Small intestine CCR6+ ILC3



Tissue Preparation: Small intestine cut into small pieces, IEL removed by washing with EDTA, dissociated with Collagenase 4, lymphocyte isolation with Percoll gradient.

Time from Mouse to Trizol: 09:15am-02:30pm


Antibody (Clone) - Fluorochrome:

CD45 (30-F11) – APC-Cy7
CD3e (145-2C11) – FITC

CD19 (eBio1D3) – FITC
Thy1.2 (53-2.1) – APC
TCR β (H57-597) – PE
CD5 (53-7.3) – PE
KLRG1 (2F1) – PE
CCR6 (140706) – BV421
NKp46 (CS96) – biotin – SA-PeCy7

ILC3: 98.4%
CCR6+: 92.1%
purity

T.8.MP.LCMV.d7.Sp: Splenic CD8 T-cell memory precursors, day 7 post LCMV



Post-sort

Tissue Preparation:

Glass slide Teasing

Time from Mouse to the 2nd sorting:

9:00am-1:00pm

Antibody (Clone) - Fluorochrome:


CCD8a(53-6.7) -- APC-Cy7

CD45.2(104) -- FITC


CD45.1(A20) -- APC

CD127(A7R34) -- PE-Cy7

KLRG1(2F1) -- Pacific Blue



T.8.Tcm.LCMV.d180.Sp: Splenic central memory CD8 T cells, day 180 post LCMV




Tissue Preparation:
Glass slide Teasing


Time from Mouse to the 2nd sorting:
9:00am-1:00pm

Antibody (Clone) - Fluorochrome:
 CCD8a(53-6.7) -- APC-Cy7
 CD45.2(104) -- FITC
 CD44(IM7) -- PE-Cy7
 CD62L(MEL-14) -- APC

Post-sort



T.8.TE.LCMV.d7.Sp: Splenic terminal effector CD8 T cells, day 8 post LCMV



Tissue Preparation:

Glass slide Teasing

Time from Mouse to the 2nd sorting:

9:00am-1:00pm

Antibody (Clone) - Fluorochrome:

CCD8a(53-6.7) -- APC-Cy7


CD45.2(104) -- FITC

CD45.1(A20) -- APC


CD127(A7R34) -- PE-Cy7

KLRG1(2F1) -- Pacific Blue

Post-sort



T.8.Tem.LCMV.d180.Sp: Splenic effector memory CD8 T cells, day 180 post LCMV



Post-sort

Tissue Preparation:

Glass slide Teasing

Time from Mouse to the 2nd sorting:

9:00am-1:00pm


Antibody (Clone) - Fluorochrome:

CCD8a(53-6.7) -- APC-Cy7


CD45.2(104) -- FITC

CD44(IM7) -- PE-Cy7


CD62L(MEL-14) -- APC



T.8.TN.P14.Sp: Splenic naive CD8 T cells from P14 TCR Tg mice



Post-sort




Tissue Preparation:
Glass slide Teasing

Time from Mouse to the 2nd sorting:
9:00am-1:00pm

Antibody (Clone) - Fluorochrome:
CD8a(53-6.7) -- FITC
CD44(IM7) -- PE-Cy7


T.8.IEL.LCMV.d7.SI: Small intestine intraepithelial lymphocytes, day7 post LCMV



Tissue Preparation:
DTE treated then isolated by Percoll gradient

Time from Mouse to the 2nd sorting:
9:00am-2:00pm

Antibody (Clone) - Fluorochrome:
CCD8a(53-6.7) -- APC-Cy7
CD45.2(104) -- FITC
CD45.1(A20) -- APC




Immunological Genome Project


B.T1.Sp: Splenic T1 B Cells

HARDY/HAYAKAWA LAB
FOX CHASE CANCER
CENTER

Pre-sort



Post-sort



Tissue Preparation: Grind between slides, ACK lysis of RBCs

Time from Mouse to Trizol: 8:30AM-11:30AM

Antibody (Clone) - Fluorochrome:

CD43(S7)-FITC; CD21/35-PE; CD3(500A-A2), Ter119, Ly6c(AL21), GR1(8C5)-Cy5PE; CD93(AA4.1)-Cy55PE; IgD(JA12.5)-Cy7PE; CD23(B3B4)-APC; CD19(1D3)-Cy55APC; CD11b(Mac1)-Cy7APC; IgM(331.12)-Cascade Blue; CD45R(6B2)-biotin/Q Dot 605 Avidin


T1 cells 99.6% purity

Immunological Genome Project


B.T2.Sp: Splenic T2 B Cells

HARDY/HAYAKAWA LAB
FOX CHASE CANCER
CENTER

Pre-sort



Post-sort



Tissue Preparation: Grind between slides, ACK lysis of RBCs

Time from Mouse to Trizol: 8:30AM-11:30AM

Antibody (Clone) - Fluorochrome:

CD43(S7)-FITC; CD21/35-PE; CD3(500A-A2), Ter119, Ly6c(AL21), GR1(8C5)-Cy5PE; CD93(AA4.1)-Cy55PE; IgD(JA12.5)-Cy7PE; CD23(B3B4)-APC; CD19(1D3)-Cy55APC; CD11b(Mac1)-Cy7APC; IgM(331.12)-Cascade Blue; CD45R(6B2)-biotin/Q Dot 605 Avidin


T2 cells 96.5% purity, T3 cells 96% purity

Immunological Genome Project


B.T3.Sp: Splenic T3 B Cells

HARDY/HAYAKAWA LAB
FOX CHASE CANCER
CENTER

Pre-sort



Post-sort



Tissue Preparation: Grind between slides, ACK lysis of RBCs

Time from Mouse to Trizol: 8:30AM-11:30AM

Antibody (Clone) - Fluorochrome:


CD43(S7)-FITC; CD21/35-PE; CD3(500A-A2), Ter119, Ly6c(AL21), GR1(8C5)-Cy5PE; CD93(AA4.1)-Cy55PE; IgD(JA12.5)-Cy7PE; CD23(B3B4)-APC; CD19(1D3)-Cy55APC; CD11b(Mac1)-Cy7APC; IgM(331.12)-Cascade Blue; CD45R(6B2)-biotin/Q Dot 605 Avidin

T2 cells 96.5% purity, T3 cells 96% purity


Immunological Genome Project

HARDY/HAYAKAWA LAB
FOX CHASE CANCER
CENTER

B1b.PC: Peritoneal Cavity B1b B Cells



Post-sort



Tissue Preparation: Peritoneal cavity washout

Time from Mouse to Trizol: 8:30AM-10:30AM


Antibody (Clone) - Fluorochrome:

CD43(S7)-FITC; CD5(53-7)-PE; CD3(500A-A2), Ter119, Ly6c(AL21), GR1(8C5)-Cy5PE; IgD(JA12.5)-Cy7PE; CD23(B3B4)-APC; CD19(1D3)-Cy55APC; CD11b(Mac1)-Cy7APC; IgM(331.12)-Cascade Blue; CD45R(6B2)-biotin/Q Dot 605 Avidin


B1b cells 93.3% purity

preB.FrD.BM: Bone marrow pre-B cells (Fraction D)

Pre-sort



Post-sort



Tissue Preparation: Harvest cells from BM using 1ml syringe, ACK lysis of RBCs


Time from Mouse to Trizol: 8:30AM-1:30PM

Antibody (Clone) - Fluorochrome:


CD43(S7)-FITC; CD127(SB/199)-PE; CD3(500A-A2), Ter119, Ly6c(AL21); CD93(AA4.1)-Cy55PE; CD24(30F1)-Alexa594; CD117(2B8)-APC; CD19(1D3)-Cy55APC; GR1(8c5), CD11b(Mac1)-Cy7APC; IgM(331.12)-Cascade Blue; CD45R(6B2)-biotin/Q Dot 605 Avidin

FrD cells 94.5% purity

proB.CLP.BM: Bone marrow pro-B cells (common lymphoid progenitors)



Post-sort



Tissue Preparation: Harvest cells from BM using 1ml syringe, ACK lysis of RBCs


Time from Mouse to Trizol: 8:30AM-1:30PM

Antibody (Clone) - Fluorochrome:


CD43(S7)-FITC; CD127(SB/199)-PE; CD3(500A-A2), Ter119, Ly6c(AL21); CD93(AA4.1)-Cy55PE; CD24(30F1)-Alexa594; CD117(2B8)-APC; CD19(1D3)-Cy55APC; GR1(8c5), CD11b(Mac1)-Cy7APC; IgM(331.12)-Cascade Blue; CD45R(6B2)-biotin/Q Dot 605 Avidin

CLP cells 91.3% purity

proB.FrA.BM: Bone marrow pro-B cells (Fraction A)



Post-sort



Tissue Preparation: Harvest cells from BM using 1ml syringe, ACK lysis of RBCs


Time from Mouse to Trizol: 8:30AM-1:30PM

Antibody (Clone) - Fluorochrome:


CD43(S7)-FITC; CD127(SB/199)-PE; CD3(500A-A2), Ter119, Ly6c(AL21); CD93(AA4.1)-Cy55PE; CD24(30F1)-Alexa594; CD117(2B8)-APC; CD19(1D3)-Cy55APC; GR1(8c5), CD11b(Mac1)-Cy7APC; IgM(331.12)-Cascade Blue; CD45R(6B2)-biotin/Q Dot 605 Avidin

FrA cells 94.3% purity

proB.FrBC.BM: Bone marrow pro-B cells (Fraction B/C)



Post-sort



Tissue Preparation: Harvest cells from BM using
1ml syringe, ACK lysis of RBCs

Time from Mouse to Trizol: 8:30AM-1:30PM

Antibody (Clone) - Fluorochrome:


CD43(S7)-FITC; CD127(SB/199)-PE; CD3(500A-A2), Ter119, Ly6c(AL21); CD93(AA4.1)-Cy55PE;
CD24(30F1)-Alexa594; CD117(2B8)-APC;
CD19(1D3)-Cy55APC; GR1(8c5), CD11b(Mac1)-
Cy7APC; IgM(331.12)-Cascade Blue;
CD45R(6B2)-biotin/Q Dot 605 Avidin

FrB/C cells 98.2% purity


Immunological Genome Project

HARDY/HAYAKAWA LAB
FOX CHASE CANCER
CENTER

B.Fo.Sp: Splenic Follicular B Cells



Post-sort



Tissue Preparation: Grind between slides, ACK lysis of RBCs

Time from Mouse to Trizol: 8:30AM-11:30AM


Antibody (Clone) - Fluorochrome:

CD43(S7)-FITC; CD21/35-PE; CD3(500A-A2), Ter119, Ly6c(AL21), GR1(8C5)-Cy5PE; CD93(AA4.1)-Cy55PE; IgD(JA12.5)-Cy7PE; CD23(B3B4)-APC; CD19(1D3)-Cy55APC; CD11b(Mac1)-Cy7APC; IgM(331.12)-Cascade Blue; CD45R(6B2)-biotin/Q Dot 605 Avidin


Follicular cells 99% purity

B.FrE.BM: Bone marrow immature B cells (Fraction E)

Pre-sort



Post-sort



Tissue Preparation: Harvest cells from BM using 1ml syringe, ACK lysis of RBCs

Time from Mouse to Trizol: 8:30AM-1:30PM

Antibody (Clone) - Fluorochrome:

CD43(S7)-FITC; CD127(SB/199)-PE; CD3(500A-A2), Ter119, Ly6c(AL21); CD93(AA4.1)-Cy55PE; CD24(30F1)-Alexa594; CD117(2B8)-APC; CD19(1D3)-Cy55APC; GR1(8c5), CD11b(Mac1)-Cy7APC; IgM(331.12)-Cascade Blue; CD45R(6B2)-biotin/Q Dot 605 Avidin


FrE cells 95.2% purity

Immunological Genome Project


HARDY/HAYAKAWA LAB
FOX CHASE CANCER
CENTER

B.MZ.Sp: Splenic Marginal Zone B Cells

Pre-sort



Post-sort



Tissue Preparation: Grind between slides, ACK lysis of RBCs

Time from Mouse to Trizol: 8:30AM-11:30AM


Antibody (Clone) - Fluorochrome:

CD43(S7)-FITC; CD21/35-PE; CD3(500A-A2), Ter119, Ly6c(AL21), GR1(8C5)-Cy5PE; CD93(AA4.1)-Cy55PE; IgD(JA12.5)-Cy7PE; CD23(B3B4)-APC; CD19(1D3)-Cy55APC; CD11b(Mac1)-Cy7APC; IgM(331.12)-Cascade Blue; CD45R(6B2)-biotin/Q Dot 605 Avidin


Marginal Zone cells 96.9% purity, B1a cells 90.9% purity

Tgd.g2+d1.24a+.Th: TCRvg2+ Scart2- immature thymocytes

First-sort



Second-sort



Tissue Preparation: Glass Tissue Grinder

Time from Mouse to Trizol: 6:00AM-4:00PM


Antibody (Clone) - Fluorochrome:

CD4 (RM4-5) PE-Cy7
 CD8(53-6.7) PE-Cy7
 CD3e (145-2C11) PerCP-Cy5.5
 TCRd (GL3) Brilliant Violet 421
 Vg2 (UC3-10A6) FITC
 Vg1.1 (2.11) APC
 CD24 (M1/69) APC-eFluor780
 Scart2 (25A2) Purified


gd.g2+d1.24a+.Th
99% purity

Tgd.g2+d1.LN: TCRvg2+ CD27+ lymphocytes

First-sort



Second-sort



Tissue Preparation: Glass Tissue Grinder

Time from Mouse to Trizol: 6:00AM-4:00PM


Antibody (Clone) - Fluorochrome:

CD4 (RM4-5) PE-Cy7
CD8(53-6.7) PE-Cy7
CD3e (145-2C11) PerCP-Cy5.5
TCRd (GL3) Brilliant Violet 421
Vg2 (UC3-10A6) FITC
Vg1.1 (2.11) APC
CD24 (M1/69) APC-eFluor780
CD27 (LG.7F9) PE


Tgd.g2+d1.LN
99% purity

Tgd.g2+d17.24a+.Th: TCRvg2+ Scart2+ immature thymocytes

First-sort



Second-sort



Tissue Preparation: Glass Tissue Grinder

Time from Mouse to Trizol: 6:00AM-4:00PM


Antibody (Clone) - Fluorochrome:

CD4 (RM4-5) PE-Cy7
CD8(53-6.7) PE-Cy7
CD3e (145-2C11) PerCP-Cy5.5
TCRd (GL3) Brilliant Violet 421
Vg2 (UC3-10A6) FITC
Vg1.1 (2.11) APC
CD24 (M1/69) APC-eFluor780
Scart2 (25A2) Purified


Tgd.g2+d17.24a+.Th
99% purity

Tgd.g2+d17.LN: TCRvg2+ CD27- lymphocytes

First-sort



Second-sort



Tissue Preparation: Glass Tissue Grinder

Time from Mouse to Trizol: 6:00AM-4:00PM


Antibody (Clone) - Fluorochrome:

CD4 (RM4-5) PE-Cy7
CD8(53-6.7) PE-Cy7
CD3e (145-2C11) PerCP-Cy5.5
TCRd (GL3) Brilliant Violet 421
Vg2 (UC3-10A6) FITC
Vg1.1 (2.11) APC
CD24 (M1/69) APC-eFluor780
CD27 (LG.7F9) PE


Tgd.g2+d17.LN
99% purity

Tgd.Sp: Splenic gdT cells

First-sort



Second-sort



Tissue Preparation: Glass Tissue Grinder

Time from Mouse to Trizol: 6:00AM-4:00PM

Antibody (Clone) - Fluorochrome:

CD4 (RM4-5) PE-Cy7

CD8(53-6.7) PE-Cy7


CD3e (145-2C11) PerCP-Cy5.5

TCRd (GL3) Brilliant Violet 421


B6.Tgd.Sp
99% purity

Tgd.g1.1+d1.24a+.Th: TCRvg1.1+ Scart2- immature thymocytes

First-sort



Second-sort



Tissue Preparation: Glass Tissue Grinder

Time from Mouse to Trizol: 6:00AM-4:00PM


Antibody (Clone) - Fluorochrome:

CD4 (RM4-5) PE-Cy7
CD8(53-6.7) PE-Cy7
CD3e (145-2C11) PerCP-Cy5.5
TCRd (GL3) Brilliant Violet 421
Vg2 (UC3-10A6) FITC
Vg1.1 (2.11) APC
CD24 (M1/69) APC-eFluor780
Scart2 (25A2) Purified


gd.g1.1+d1.24a+.Th
99% purity

Tgd.g1.1+d1.LN: TCRvg1.1+ CD27+ lymphocytes

First-sort



Second-sort



Tissue Preparation: Glass Tissue Grinder


Time from Mouse to Trizol: 6:00AM-4:00PM

Antibody (Clone) - Fluorochrome:

CD4 (RM4-5) PE-Cy7
CD8(53-6.7) PE-Cy7
CD3e (145-2C11) PerCP-Cy5.5
TCRd (GL3) Brilliant Violet 421
Vg2 (UC3-10A6) FITC
Vg1.1 (2.11) APC
CD24 (M1/69) APC-eFluor780
CD27 (LG.7F9) PE

Tgd.g1.1+d1.LN
99% purity

NK.27+11b+.Sp: Splenic semi-mature NK cells



Tissue Preparation: [Splenocytes](#)


Time from Mouse to TCL buffer (1% 2-ME) or BAMBANKER:

7:00 am – 2-3:00 pm

Antibody (Clone) - Fluorochrome:

CD49b (DX5) – FITC
NK1.1 (PK136) – PerCP-Cy5.5
CD3e (145-2C11) – PE-Cy7
CD19 (1D3) – PE-Cy7
CD11b (M1/70) – PE
CD27 (LG.3A10) – APC
CD49a (REA493) – Biotin
CD127 (A7R34) – Biotin
CD51 (RMV-7) – Biotin
Streptavidin-BV421

NK.27-11b+.BM: Bone marrow mature NK cells




Tissue Preparation: [Splenocytes](#)

Time from Mouse to TCL buffer (1% 2-ME) or BAMBANKER:


7:00 am – 2-3:00 pm

Antibody (Clone) - Fluorochrome:

CD49b (DX5) – FITC
NK1.1 (PK136) – PerCP-Cy5.5
CD3e (145-2C11) – PE-Cy7
CD19 (1D3) – PE-Cy7
CD11b (M1/70) – PE
CD27 (LG.3A10) – APC
CD49a (REA493) – Biotin
CD127 (A7R34) – Biotin
CD51 (RMV-7) – Biotin
Streptavidin-BV421



NK.27-11b+.Sp: Splenic mature NK cells




Tissue Preparation: [Splenocytes](#)

Time from Mouse to TCL buffer (1% 2-ME) or BAMBANKER:


7:00 am – 2-3:00 pm

Antibody (Clone) - Fluorochrome:

CD49b (DX5) – FITC
NK1.1 (PK136) – PerCP-Cy5.5
CD3e (145-2C11) – PE-Cy7
CD19 (1D3) – PE-Cy7
CD11b (M1/70) – PE
CD27 (LG.3A10) – APC
CD49a (REA493) – Biotin
CD127 (A7R34) – Biotin
CD51 (RMV-7) – Biotin
Streptavidin-BV421



NK.27+11b-.BM: Bone marrow immature NK cells




Tissue Preparation: [Splenocytes](#)

Time from Mouse to TCL buffer (1% 2-ME) or BAMBANKER:


7:00 am – 2-3:00 pm

Antibody (Clone) - Fluorochrome:

CD49b (DX5) – FITC
NK1.1 (PK136) – PerCP-Cy5.5
CD3e (145-2C11) – PE-Cy7
CD19 (1D3) – PE-Cy7
CD11b (M1/70) – PE
CD27 (LG.3A10) – APC
CD49a (REA493) – Biotin
CD127 (A7R34) – Biotin
CD51 (RMV-7) – Biotin
Streptavidin-BV421



NK.27+11b-.Sp: Splenic immature NK cells




Tissue Preparation: [Splenocytes](#)

Time from Mouse to TCL buffer (1% 2-ME) or BAMBANKER:


7:00 am – 2-3:00 pm

Antibody (Clone) - Fluorochrome:

CD49b (DX5) – FITC
NK1.1 (PK136) – PerCP-Cy5.5
CD3e (145-2C11) – PE-Cy7
CD19 (1D3) – PE-Cy7
CD11b (M1/70) – PE
CD27 (LG.3A10) – APC
CD49a (REA493) – Biotin
CD127 (A7R34) – Biotin
CD51 (RMV-7) – Biotin
Streptavidin-BV421



NK.27+11b+.BM: Bone marrow semi-mature NK cells



Tissue Preparation: [Splenocytes](#)


Time from Mouse to TCL buffer (1% 2-ME) or BAMBANKER:

7:00 am – 2-3:00 pm

Antibody (Clone) - Fluorochrome:


CD49b (DX5) – FITC
NK1.1 (PK136) – PerCP-Cy5.5
CD3e (145-2C11) – PE-Cy7
CD19 (1D3) – PE-Cy7
CD11b (M1/70) – PE
CD27 (LG.3A10) – APC
CD49a (REA493) – Biotin
CD127 (A7R34) – Biotin
CD51 (RMV-7) – Biotin
Streptavidin-BV421

Post-sort




DC.8+.Sp: Splenic CD8+ Dendritic cells

Pre-sort



Post-sort



Tissue Preparation: Collagenase V digestion, push through 70um cell strainer


Time from Mouse to Trizol: 09:00am-12:00pm

Antibody (Clone) - Fluorochrome:

- CD3, CD19 – FITC
- CD8 – PE
- CD4 – APC
- CD11c – PE/Cy7
- CD45 – BV510
- MHC II – APC/Cy7
- viability - PI

10000 cells
98.0% purity

DC.pDC.Sp: Splenic plasmacytoid dendritic cells




Tissue Preparation: Filter/syringe plunger teasing

Time from Mouse to Trizol: 09:30am-01:00pm


Antibody (Clone) - Fluorochrome:

CD3 (eBio500A2)
CD19 (eBio1D3) FITC
CD11c (N418) PE
B220 (RA3-6B2) eFluor450
PDCA1 (eBio927) APC
Siglec H (eBio440c) PE-Cy7




MF.LP.SI: Small intestine lamina propria macrophages

Pre-sort



Second Sort



Tissue Preparation: Collagenase DNase


Time from Mouse to Trizol: 9:00 am-2:00pm

Antibody (Clone) - Fluorochrome:


CD45 -- APC-Cy7
CD11c -- PE-Cy7
MHCII -- PacificBlue
CD64 -- APC
CD103 -- PE
CD11b -- FITC

MF.Alv.Lu: Lung alveolar macrophages

Pre-sort



Post-sort



Tissue Preparation: Collagenase V digestion,
dounce glass homogenizer

Time from Mouse to Trizol: 09:00am-12:00am


Antibody (Clone) - Fluorochrome:

Siglec F – PE
CD11c (N418) – PE/Cy7
CD45 (30-F11) – BV510
viability - PI


10000 cells
97.5% purity

MF.plC.Alv.Lu: Poly(I:C)-stimulated lung alveolar macrophages

Pre-sort



Post-sort



Tissue Preparation: Collagenase V digestion,
dounce glass homogenizer

Time from Mouse to Trizol: 09:00am-12:00am


Antibody (Clone) - Fluorochrome:

Siglec F – PE
CD11c – PE/Cy7
CD45 – BV510
viability - PI


10000 cells
96.5% purity

DC.103+11b-.SI: Small intestine CD103+ CD11b- dendritic cells

Pre-sort



Second Sort



Tissue Preparation: [Collagenase DNase](#)


Time from Mouse to Trizol: [9:00 am-2:00pm](#)

Antibody (Clone) - Fluorochrome:

CD45 -- APC-Cy7
CD11c -- PE-Cy7
MHCI -- PacificBlue
CD64 -- APC
CD103 -- PE
CD11b -- FITC

DC.103+11b+.SI: Small intestine CD103+ CD11b+ dendritic cells

Pre-sort




Second Sort

Tissue Preparation: Collagenase DNase

Time from Mouse to Trizol: 9:00 am-2:00pm


Antibody (Clone) - Fluorochrome:

CD45 -- APC-Cy7
CD11c -- PE-Cy7
MHCII -- PacificBlue
CD64 -- APC
CD103 -- PE
CD11b -- FITC




DC.4+.Sp: Splenic CD4+ Dendritic cells

Pre-sort



Post-sort



Tissue Preparation: Collagenase V digestion, push through 70um cell strainer

Time from Mouse to Trizol: 09:00am-12:00pm

Antibody (Clone) - Fluorochrome:

CD3, CD19 - FITC

CD4 - APC

CD11c - PE/Cy7

CD8 - PE

CD45 - BV510


MHC II - APC/Cy7

viability - PI


10000 cells
97.1% purity

GN.BM: Bone marrow neutrophils

Pre-sort



2nd sorting



Tissue Preparation:

Flushed BM from femur using syringe, ACK lysis

Time from Mouse to the 2nd sorting:


8am to 11am

Antibody (Clone) - Fluorochrome:


CD11b (M1/70) – PE

Ly6-G (RB6-8C5) – APC

GN.SP: Splenic neutrophils



Post-sort




Tissue Preparation: Mechanical disruption of spleen
Time from Mouse to TCL Buffer: 8:30am-11:00am


Antibody (Clone) - Fluorochrome:

TCRb (H57-597)—Pacific Blue
B220 (RA3-6B2)—Pacific Blue
CD11c (N418)—Alexa 700
MHCII (M5/114.15.2)—FITC
Ly6G (RB6-8C5)—APC

GN.Thio.PC: Thio-induced peritoneal cavity neutrophils



2nd sorting



Tissue Preparation:

1ml 3% Thio(SIGMA, in dw) ip injected 18hrs before sac.
Peritoneal lavage with Phenol-red free DMEM, ACK lysis

Time from Mouse to the 2nd sorting:


8am to 10am

Antibody (Clone) - Fluorochrome:


CD11b (M1/70) – PE

Ly6-G (RB6-8C5) – APC-Cy7

B.PB.Sp: Splenic plasmablasts



Second sort



Tissue Preparation: Bones crushed with mortar and pestle, ACK lysis of RBC
Time from Mouse to Trizol: 4 hours

Antibody (Clone) - Fluorochrome:

CD138 (281.2) – PE-Cy7

NK1.1 (PK136) – APC


CD11b (M1/70) – APC

TCRb (H57-597) – APC


MHCII () – APC-Cy7

Over 90 %
before resort

B.PC.BM: Bone marrow plasma cells



Second sort




Tissue Preparation: Bones crushed with mortar and pestle, ACK lysis of RBC
Time from Mouse to Trizol: 4 hours

Antibody (Clone) - Fluorochrome:


CD138 (281.2) – PE-Cy7
NK1.1 (PK136) – APC
CD11b (M1/70) – APC
TCRb (H57-597) – APC

Over 90 %
before resort

B.PC.Sp: Splenic plasma cells



Second sort




Tissue Preparation: Spleen homogenized through 100 μ m cell-strainer, ACK lysis
Time from Mouse to Trizol: 4 hours

Antibody (Clone) - Fluorochrome:

CD138 (281.2) – PE-Cy7
NK1.1 (PK136) – APC
CD11b (M1/70) – APC
TCRb (H57-597) – APC
MHCII () – APC-Cy7

Over 95 %
before resort

B.GC.CB.Sp: Splenic germinal center centroblasts




Second sort

Tissue Preparation: Spleen homogenized through 100 μ m cell-strainer, ACK lysis, positive selection with anti-IgG1-microbeads (Miltenyi)

Time from Mouse to Trizol: 5 hours


Antibody (Clone) - Fluorochrome:

CD138 (281.2) – PE
Gr1 (RB6-6C5) – PE
TCRb (H57-597) – PE
B220 (RA3-6B2) – FITC
Fas (Jo2) – PE-Cy7
CD38 (90) – AF780
CXCR4 (2B11) – BV421
CD86 (GL1) – BV605




Over 90 %
before resort

B.GC.CC.Sp: Splenic germinal center centrocytes



Second sort



Tissue Preparation: Spleen homogenized through 100 μ m cell-strainer, ACK lysis, positive selection with anti-IgG1-microbeads (Miltenyi)


Time from Mouse to Trizol: 5 hours

Antibody (Clone) - Fluorochrome:


CD138 (281.2) – PE
Gr1 (RB6-6C5) – PE
TCRb (H57-597) – PE
B220 (RA3-6B2) – FITC
Fas (Jo2) – PE-Cy7
CD38 (90) – AF780
CXCR4 (2B11) – BV421
CD86 (GL1) – BV605

Over 90 %
before resort

B.mem.Sp: Splenic memory B cells



Second sort



Tissue Preparation: Spleen homogenized through 100 μ m cell-strainer, ACK lysis, positive selection with anti-IgG1-microbeads (Miltenyi)



Time from Mouse to Trizol: 5 hours

Antibody (Clone) - Fluorochrome:

CD138 (281.2) – PE
Gr1 (RB6-6C5) – PE
TCRb (H57-597) – PE
B220 (RA3-6B2) – FITC
Fas (Jo2) – PE-Cy7
CD38 (90) – AF780
IgM (II/41) – APC-eF780
IgD (11-26C) – APC-eF780
IgG (poly) – Cy5

Over 60 %
before resort

MF.ICAM+480hi.PC: Large peritoneal cavity macrophages




Tissue Preparation:
Peritoneal lavage

Time from Mouse
30min antibody staining

Antibody (Clone) - Fluorochrome
 ICAM2 -- Alexa 488
 CD11b -- APC-Cy7
 CD115 -- APC
 MHCII -- efluor450
 F4/80 -- PE-Cy7

10,000 cells
Purity >95%


MF.microglia.CNS: Microglias



Tissue Preparation:
Percoll


Time from Mouse
20min antibody staining

Antibody (Clone) - Fluorochrome
 CD19, CD3, Ly6G -- PE
 CD11b -- APC-Cy7
 CD64 -- APC
 CD45 -- efluor450
 CD206 -- FITC
 F4/80 -- PE-Cy7



10,000 cells
Purity >95%


MF.RP.Sp: Spleen red pulp macrophages



Tissue Preparation:
Mechanical dissociation



Time from Mouse
30min antibody staining

Antibody (Clone) - Fluorochrome
CD19, Ly6G, CD115 -- PE
CD11b -- APC-Cy7
CD64 -- APC
CD45 -- efluor450
MerTK -- PE-Cy7
F4/80 -- Alexa488



10,000 cells
Purity >95%

Mo.6C+II-.Bl: Ly6Chi blood monocytes





Tissue Preparation:
500uL of Lysis buffer

Time from Mouse
20min antibody staining

Antibody (Clone) - Fluorochrome:
B220, CD3, Ly6G -- PE-Cy7
CD11b -- APC-Cy7
CD115 -- APC
Ly6C -- efluor450

10,000 cells
Purity>95%

Mo.6C-II-.Bl: Ly6Clo blood monocytes




Tissue Preparation:
500uL of Lysis buffer

Time from Mouse
20min antibody staining

Antibody (Clone) - Fluorochrome
B220, CD3, Ly6G -- PE-Cy7
CD11b -- APC-Cy7
CD115 -- APC
Ly6C -- efluor450

10,000 cells
Purity>95%

MF.226+II+480lo.PC: Small peritoneal cavity macrophages



Tissue Preparation:
Peritoneal lavage

Time from Mouse
20min antibody staining

Antibody (Clone) - Fluorochrome

ICAM2 -- Alexa 488

CD11b -- APC-Cy7


CD115 -- APC

MHCII -- efluor450


CD226 -- PE

10,000 cells
Purity >95%

B.Fem.SP: Splenic B cell (Female)



Post-sort



Tissue Preparation: Mechanical disruption of spleen

Time from Mouse to TCL Buffer: 10:30am-12:00pm


Antibody (Clone) - Fluorochrome:

TCRb (H57-597)—Pacific Blue


CD19 (1D3)—PE-Cy7

IgM (eB121-15F9)—FITC

MF.Fem.PC: Peritoneal cavity macrophages (female)



2nd sorting



Tissue Preparation:

Peritoneal lavage with Phenol-red free DMEM, ACK lysis

Time from Mouse to the 2nd sorting:

8am to 10am


Antibody (Clone) - Fluorochrome:

F4/80 (BM8) -- PE-Cy7


ICAM2 (3C4 (mIC2/4)) -- FITC

T.4.Nve.Fem.Sp: Splenic Naïve CD4 T cells (female)

Pre-sort



2nd sorting



Tissue Preparation:

100um mesh/syringe plunger teasing, ACK lysis


Time from Mouse to the 2nd sorting:

8am to 10am

Antibody (Clone) - Fluorochrome:

CD4 (RM4-5) – PE, CD8a (53-6.7) – A700, TCRb (H57-597) – eF450,
CD44(1M7) – FITC, CD62L(Mel-14) – PE-Cy7, CD25(PC61) -- APC
dump = CD19(1D3), Ter119(Ter-119) , Gr1(RB6-8C5) – APC-eF780
NK1.1(PK136), F4/80(BM8) -- APC-Cy7

LEC.SLN: Lymphatic Endothelial cells



Tissue Preparation:


Dispase and Collagenase P digestion with DNaseI, followed by a depletion employing CD45 microbeads (MACS)

Time from Mouse to the 2nd sorting:

9am to 4pm


Antibody (Clone) - Fluorochrome:

PDPN (8.1.1) -- APC
CD31 (390) -- PE/DAZZLE 594
CD45 (30-F11), Ter119 (TER119), Epcam (G8.8) -- PE-Cy7
Calcein Blue
Live/Dead -- Near-IR633



Purity is
90.1% of
cell gate

BEC.SLN: Blood Endothelial cells



Tissue Preparation:


Dispase and Collagenase P digestion with DNaseI, followed by a depletion employing CD45 microbeads (MACS)

Time from Mouse to the 2nd sorting:

9am to 4pm


Antibody (Clone) - Fluorochrome:

PDPN (8.1.1) -- APC
CD31 (390) -- PE/DAZZLE 594
CD45 (30-F11), Ter119 (TER119), Epcam (G8.8) -- PE-Cy7
Calcein Blue
Live/Dead -- Near-IR633



Purity is
91.4% of
cell gate

FRC.SLN: Fibroblastic reticular cells



Tissue Preparation:


Dispase and Collagenase P digestion with DNaseI, followed by a depletion employing CD45 microbeads (MACS)

Time from Mouse to the 2nd sorting:

9am to 4pm


Antibody (Clone) - Fluorochrome:

madcam1 (MECA-367) -- Alexa488
CD21/35 (7E9) -- FITC
PDPN (8.1.1) -- APC
CD31 (390) -- PE/DAZZLE 594
CD45 (30-F11), Ter119 (TER119), Epcam (G8.8) -- PE-Cy7
Calcein Blue
Live/Dead -- Near-IR633



Purity is
98.5% of
cell gate

IAP.SLN: Pericytes



Tissue Preparation:


Dispase and Collagenase P digestion with DNaseI, followed by a depletion employing CD45 microbeads (MACS)

Time from Mouse to the 2nd sorting:

9am to 4pm

Antibody (Clone) - Fluorochrome:

PDPN (8.1.1) -- APC
CD31 (390) -- PE/DAZZLE 594
CD45 (30-F11), Ter119 (TER119), Epcam (G8.8) -- PE-Cy7
Calcein Blue
Live/Dead -- Near-IR633



Purity is
99.2% of
cell gate