




Tissue preparation: Glass slide teasing, ACK lysis

Time from mouse to Trizol: 8.45am-1pm

Antibody (Clone) - Fluorescence:

CD11b (M1/70, CD11c (HL3), CD19 (MB19.1), CD49b (HMA2), Gr1 (RB6-8C5), Ter119 (TER-119) - Biotin - SAPEAlexa610

CD4 (RM4-5) - APC

CD8a (CT-CD8a) - FITC

CD69 (H1.2F3) - PE

TCRb (H57-597) - APCCy7